Questionnaire – FACS Sorting (BD FACSAria IIu)

| Experimentalist: | | | |
|-------------------------|-------------------|----------------|----------------------|
| Principal Investigator: | | | |
| Institute: | | | |
| Phone: | | | |
| e-mail: | | | |
| Sorting Material: | | | |
| Species of origin: | □human | | |
| | ∐humaı | n cell line | |
| | ∏primaı | y tissue | |
| | □mouse | □other | |
| Organ of origin: | ПРВ | □spleen | |
| | □LN | · | |
| Condition: | □non-fixed | ∏fixed | |
| | □sterile | □non-sterile |) |
| | □virally trans | duced or infed | eted |
| | specify | | |
| Fluorescent Dyes: | □FITC | □PE | □CY3 |
| | □PerCP | □CY5 | □APC |
| | □PacificBlue | □7AAD | □PI |
| | □GFP | | |
| | □other | | |
| Sort Parameters | □normal | | ☐single cell deposit |
| | tubes | | tubes |
| | well pla | ate | well plate |
| Temperature Control | ☐room temperature | | □4°C |

Important information for cell sorting at the FACSAria Ilu

- To schedule your sort appointments call 5549 (Mr. Eder or Ms. Fink) at least **ONE WEEK** prior to intended sort date
- Appointments can be booked Monday to Friday between 8 am and 4 pm (other times upon request)
- Please provide a flow cytometric analysis of the intended sort material (or one as similar as possible) and discuss it either with one of the operators or with PD Dr. Petra Hoffmann (phone: 5592) prior to the first sort date
- If you are late for your appointment please give us a call (lab: 5549 or at the sorter: 6644). Delay may result in reduced time for sorting
- Please provide a telephone number for contact on the sort day
- Ideal cell concentration for most sorting material: 20x10⁶ cells/ml
- Minimum sample volume: 0,5 ml
- Sorting material must be prepared and kept under sterile conditions (exceptions only upon prior request)
- Cells must be resuspended in sterile PBS (addition of serum/protein possible: please consult operators) and filtered immediately before sorting (≤ 40µm)
- **PLEASE NOTE**: Safety level **S1** material (genetic engineering and biohazard) can be sorted without prior notice. Upon special request and ONLY AFTER PRIOR CONSULTATION we might be able to sort material with biohazard level S2.

| no infectious or otherwise hazardous materials for man. | | |
|---|------------------------|--|
| | | |
| | | |
| Date | Principal Investigator | |

I herewith certify that the cells provided by me or my member of staff contain